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**Notes on the interpretation of EDX spectra and result tables:****Englisch**

Energy-dispersive microanalysis (short: EDX) is a method for surface-analysis. In addition to the advantage of being able to carry out analyses locally and from small surface areas, there are the following disadvantages that must be taken into account when observing the spectra and interpreting the results:

1. The method is only suitable for the **detection of elements**, not compounds.
2. **Line overlays** (Certain elements cannot be distinguished from each other because their line positions are too close together. Based on his experience, the device operator has usually decided in each case for the element that is likely to occur in the spectrum designation, if line overlays occurred. However, it cannot be excluded that one or more of the other possible elements are also present in very low concentrations).
3. **Detection limit** (Depending on the atomic number of an element, the excitation voltage during the examination and other factors, the detection of low element concentrations, approx. <0.1wt%, is not possible. In particular, light elements such as carbon or oxygen are difficult to detect quantitatively because line layers of other elements in the low-energy range often overlay these element lines.)
4. **Excitation conditions** (Smallest particles can be "shot through" by the excitation beam, so that the analysis results may also originate from the particle neighbourhood or the subsurface. Furthermore, in the case of layered structures, individual thin layers are usually indistinguishable from each other because the excitation takes place in a volume spanning all layers).
5. **Quantification** (A standard-free quantification is carried out according to a basic physical model with the ESPRIT software from Bruker, whereby elements such as C, O can be determined in difference to 100%. Standard-free quantification requires a homogeneous distribution of the elements present in the excited sample volume. The quantification of rough and heterogeneous surfaces is therefore highly error-prone! Oxide layers on surfaces also contribute to the heterogeneity of the excited sample volume. The result of a standard-free quantitative analysis represents an estimate in the sense of ISO22309).

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